



hello@belloeveryday.com



**bello**  
Digital Belly Fat Scanner

# Quick-Start Guide

Revision 2.0

© 2020 Olive Healthcare Inc. All rights reserved.

## Bello Mobile App Installation

### iOS Users

- 1 Enter the App Store URL  
<https://apps.apple.com/app/id1485202681>
- 2 Search App Stores by Keyword  
#Bello #Bello Decode

- 3 Scan QR Code



### Android Users

- 1 Enter the App Store URL  
<https://play.google.com/store/apps/details?id=ohc.app.bello>
- 2 Search App Stores by Keyword  
#Bello #Bello Decode

- 3 Scan QR Code



## Bello Bluetooth Pairing


- 1 Sign up or log in on the Bello app. Once you enter your basic information and get to the "Connect your Bello" screen, **turn on your device** (ref. #03) and wait until it shows up on the app screen.
- 3 You are ready for your first scan when the LED checkmark appears on your Bello device and the app shows the "Scan upper belly" screen

- 2 All iOS/Android users have to click on the yellow "Connect" button and enter **6-digit** PIN Code "111111" on the pop-up screen.

**▲ Caution :** All users have to enter **6-digit** PIN Code "111111" as stated above. Android users will not be able to pair with their Bello device if they enter a 4-digit code such as "0000", "1111", or "1234". If you have multiple Bello devices turned on near you, the device that is closest to you is usually shown at the top of the device list.


## Turning the device ON or OFF

### 1 Turning the device ON



Press the top button 2 seconds to turn on the device. You will hear a **"beep"** sound and device will stay on for approximately 3 minutes. The device will turn off automatically if not used.


### 2 Turning the device OFF




To turn off the device, press the same button for more than 3 seconds.

**⚠ Caution :** If you press and hold the top button down for 10 seconds for a factory reset, all bluetooth pairing data on the Bello device will be lost and users will have to repeat the initial bluetooth pairing process to connect to the Bello.

## How to scan your belly fat



When the "Scan upper belly" screen appears on the app, place the Bello device one inch above your belly button and **lightly tap the top button**.



The measurement will take approximately 3 to 5 seconds and you will hear a **"beep"**, indicating the upper belly scan is complete.

Continue with "Scan lower belly" in the same manner.

When the lower belly scan is completed, you will receive your scan result.

### ▲ Caution :

Please hit the top button shortly. If you press and hold the top button for 2 to 3 seconds, the device is powered off.

## FCC Statement

### FCC Part 15.19

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

### FCC Part 15.21

Any changes or modifications (including the antennas) to this device that are not expressly approved by the manufacturer may void the user's authority to operate the equipment.

#### FCC RF Radiation Exposure Statement:

This equipment complies with FCC RF Radiation exposure limits set forth for an uncontrolled environment.

This device and its antenna must not be co-located or operating in conjunction with any other antenna or transmitter.

### FCC Part 15.105

**Note:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Modifications not expressly approved by the manufacturer could void your authority to operate the equipment under FCC rules.


## Bello Inquiry

### Two ways to send Bello Inquiry to get customer support

#### 1. Tap the "Email Us" button on the Profile page

- When you tap the "Email Us" button (✉️) at the top right of the screen, you will see Bello Inquiry page.
- By going on the "Bello Inquiry" page, you can send inquiry to **hello@belloeveryday.com** including some information as stated in "Caution Box" below.

#### 2. Send Bello Inquiry directly to *hello@belloeveryday.com*



iOS User



Android User

#### ▲ Caution :

Bello inquiry should include some information as stated below:

- Bello account (Email)
- Phone model
- OS version
- Bello App. version
- Your inquiry

This information will help us assist you faster.

FAQ



User Manual

